


## CHIP-SEQ DATA ANALYSIS WORKFLOW



- **Goal:**Chromatin immunoprecipitation (ChIP) experiments are performed to identify DNA that binds specific (chromatin) proteins of interest.

### Data Analysis Pipeline Steps:

- Sequence Reads
- Quality Control (Evaluate the quality of the sequencing data)à Fastqc
- Alignment to Genome & Result Output to SAM à Bowtie2
- SAM to BAM conversion à Samtools
- Sorting BAM by Genomic Coordinates and filtering only uniquely mapped reads à Sambamba
- Peak Calling (Identify **areas** in the genome that have been enriched with **aligned reads** as a consequence of performing a ChIP-Seq experiment) à MACS
- Downstream Analysis (annotation of the peaks)

**Organism:** Oryza sativa

Secondary Study Accession: SRP017137

Study Title:

Chromatin immunoprecipitation of MADS29 to find the cis regulatory regions and targets of this transcription factor

## Source

- [DNA](#)(3)
  - **Library Layout**
- [single](#)(3)
  - **Platform**
- [Illumina](#)(3)
  - **Strategy**
- [other](#)(3)

**TOTAL biosample accession no. =**

**SRR616062**

**SRR616063**

**SRR616064**

## **Project: PRJNA179281**

A ChIP-seq analysis revealed 378 targets of OsMADS29, which include genes involved in cytokinin metabolism and auxin signaling, carbohydrate metabolism, transporters and a large number of transcription factors, reflecting on its functional diversity. Overall design: Chromatin Immunoprecipitation of MADS29 using antiMADS29 antibody was done using rice 9 DAP seed tissue with total input chromatin and mock immunoprecipitated chromatin as controls

## **#FastQC**

```
sudo apt-get install fastqc
```

## **#Trimmomatic**

```
sudo apt-get update
```

```
sudo apt-get install trimmomatic-0.36
```

```
for single end=
```

```
~/Documents/chip_seq/ena_files(path)
```

```
java -jar trimmomatic-0.36.jar SE SRR616062/SRR616062.fastq.gz forward62.fq
```

```
ILLUMINACLIP:TruSeq3-SE.fa:2:30:10 LEADING:3 TRAILING:3 SLIDINGWINDOW:4:30
```

```
MINLEN:20
```

```
(output= TrimmomaticSE: Started with arguments:
```

```
SRR616062/SRR616062.fastq.gz forward62.fq ILLUMINACLIP:TruSeq3-SE.fa:2:30:10
```

```
LEADING:3 TRAILING:3 SLIDINGWINDOW:4:30 MINLEN:20
```

```
Automatically using 1 threads
```

Using Long Clipping Sequence: 'AGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGT'  
Using Long Clipping Sequence: 'AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC'  
ILLUMINACLIP: Using 0 prefix pairs, 2 forward/reverse sequences, 0 forward only sequences, 0 reverse only sequences  
Quality encoding detected as phred33  
Input Reads: 11791914 Surviving: 10076007 (85.45%) Dropped: 1715907 (14.55%)  
TrimmomaticSE: Completed successfully)  
(create a forward62.fa)

## #bowtie2

sudo apt-get install bowtie2  
The bowtie2-build indexer-  
Command Line

Usage:

#bowtie2-build [options]\* <reference\_in> <bt2\_base>  
bowtie2-build -f GCF\_001433935.1\_IRGSP-1.0\_genomic.fna.gz b2\_base  
bowtie2-build builds a Bowtie index from a set of DNA sequences. bowtie2-build outputs a set of 6 files with suffixes.  
(create file =(1) b2\_base.1.bt2 , (2) b2\_base.2.bt2, (3) b2\_base.3.bt2, (4)b2\_base.4.bt2, (5)b2\_base.rev.1.bt2 , (6)b2\_base.rev.2.bt2

#Getting started with Bowtie 2=  
bowtie2 -x b2\_base -U forward62.fq -S forward62.sam(open terminal where all index file)  
(output- 0076007 reads; of these:  
10076007 (100.00%) were unpaired; of these:  
2921040 (28.99%) aligned 0 times  
4355722 (43.23%) aligned exactly 1 time  
2799245 (27.78%) aligned >1 times

71.01% overall alignment rate)  
head forward62.sam(commands)

(outputs=  
@HD VN:1.0 SO:unsorted  
@SQ SN:NC\_029256.1 LN:43270923  
@SQ SN:NC\_029257.1 LN:35937250  
@SQ SN:NC\_029258.1 LN:36413819  
@SQ SN:NC\_029259.1 LN:35502694  
@SQ SN:NC\_029260.1 LN:29958434  
@SQ SN:NC\_029261.1 LN:31248787  
@SQ SN:NC\_029262.1 LN:29697621  
@SQ SN:NC\_029263.1 LN:28443022  
@SQ SN:NC\_029264.1 LN:23012720)  
bowtie2 -x b2\_base -U forward63.fq -S forward63.sam

(output=33812514 reads; of these:  
33812514 (100.00%) were unpaired; of these:  
2809217 (8.31%) aligned 0 times  
16956308 (50.15%) aligned exactly 1 time  
14046989 (41.54%) aligned >1 times

91.69% overall alignment rate)  
bowtie2 -x b2\_base -U forward64.fq -S forward64.sam

```
(output= e2 -x b2_base -U forward64.fq -S forward64.sam
137885 reads; of these:
  137885 (100.00%) were unpaired; of these:
    11074 (8.03%) aligned 0 times
    69871 (50.67%) aligned exactly 1 time
    56940 (41.30%) aligned >1 times
91.97% overall alignment rate)
head forward62.sam
head -20 forward62.sam(open sam file )
#SAM TO BAM
IN SAM file(SRR616062.10525368 0      NC_029256.1 5838198      42      36M      *      00
      GGAAGCTCTGGTGGAGCTGCAGCGCGTGC GCGCCGT
      IIIIIIIIIIIIIIIIGIIIIIIIIIIIIIG  AS:i:0 X)
```

### #install samtools

```
The samtools help= (samtools view ,samtools sort ,samtools depth|)
#Converting SAM to BAM with samtools "view"
samtools view -S -b sample.sam > sample.bam
samtools view sample.bam | head
#samtools "sort"
samtools view sample.bam | head
samtools sort sample.bam -o sample.sorted.bam
#Now let's check the order:
samtools view sample.sorted.bam | head
samtools view sample.sorted.bam | wc -l(#count total alignment= 10076007)
samtools view -H sample.sorted.bam(inspect header)
samtools view forward62.sorted.bam | awk '{print $5}' | sort -k1g| uniq -c
samtools view forward62.sorted.bam | awk '{print $2}' | sort | uniq -c
samtools view -f 0*2 forward62.sorted.bam |head(capture the FLAG)
samtools view -f 0*2 forward62.sorted.bam |wc -l(how many properly alignment there)
```

### #install sambamba tool(filtering bam file(remove time alignment sequence, remove >1 time alignment sequence we should only one time sequence alignment)

```
sambamba-0.8.0-linux-amd64-static.gzsambamba-0.8.0-linux-amd64-static.gz
chmod u+x sambamba-0.8.0
./sambamba-0.8.0
Available commands:
```

```
view      view contents and convert from one format
          to another (SAM/BAM/JSON/UNPACK)
index     build index (BAI)
merge     merge files (BAM)
sort      sort file (BAM)
slice     slice file (BAM using BED)
markdup   mark or remove duplicates (BAM)
subsample subsample (BAM)
flagstat  output statistics (BAM)
depth     output statistics (BAM)
validate  simple validator (BAM)
```

```
./sambamba-0.8.0 -h
./sambamba-0.8.0 view forward64.bam | head |wc -l(# total world 8868193)

./sambamba-0.8.0 view forward64.bam | head |wc
```

```
./sambamba-0.8.0 view forward64.bam | head
sambamba markdup OPTIONS <input.bam> <output.bam> (./sambamba-0.8.0 markdup -r
forward62.bam output62.bam)
./sambamba-0.8.0 markdup -r forward63.bam output63.bam(remove duplicate aligned sequence)
./sambamba-0.8.0 markdup -r forward62.bam markdup62.bam(markdup.bam output file)
(output= inding positions of the duplicate reads in the file...
sorted 0 end pairs
and 31003297 single ends (among them 0 unmatched pairs)
collecting indices of duplicate reads... done in 3580 ms
found 9919438 duplicates
collected list of positions in 0 min 32 sec
removing duplicates...
collected list of positions in 0 min 50 sec)
#MACS2(INSTALL MACS2 tools cammands and requirement tools)
Requirements
1.1.1 Software
• Python 2.7.10+
• Cython 0.25+
• HTSlib
```

### • SAMtools

```
(sudo apt-get install make build-essential zlib1g-dev libbz2-dev liblzma-dev curl pigz
, → libncurses5-dev)
pip install Cython
git clone https://github.com/samtools/htslib.git
cd htslib
autoheader
autoconf
```

### MACS2 installation

```
git clone https://github.com/Multiscale-Genomics/mg-process-macs2.git
pip install git+https://github.com/Multiscale-Genomics/mg-process-macs2.git
Usage: macs2 <-t tfile> [-n name] [-g genomesize] [options]
```

```
Example: macs2 -t ChIP.bam -c Control.bam -f BAM -g hs -n test -B -q 0.01
macs2 callpeak -t output62.bam -f BAM -g 4.3e+6 --nomodel -n test -B --outdir macs2_1 -q 0.01
samtools sort output62.bam -o output62.sorted.bam
macs2 callpeak -t uniq62.bam -f BAM --name=peak1q -q 0.01( file create) or macs2 callpeak -t
uniq62.bam -f BAM --name=peak1q --outdir macs2_62(file name) -q 0.01
Rscript NAME_model.r
Rscript peak1_model.r(create model.r)
Rscript peak1q_model.r( all output files start peak1q_)
macs2 callpeak -t uniq64.bam -f BAM --name=peak1q --outdir macs2_64 -q 0.01
Rscript peak1q_model.r
macs2 callpeak -t uniq63.bam -f BAM --name=peak1q --outdir macs2_63 -q 0.01
Rscript peak1q_model.r(open macs2_63 file then use cammand)
```